FotonaSmooth

REVOLUTIONARY LASER FOR VAGINAL HEALTH



FotonaSmooth®: In-Office Based Non-Invasive Laser Therapy for Vaginal Health

Pelvic floor dysfunction

Symptoms of age- and childbirth-related pelvic floor dysfunction, such as urinary incontinence and pelvic floor prolapse, affect hundreds of millions of women worldwide.

- 60-80% of women over 50 will experience atrophy in their lifetime
- 40% of women suffer from some form of urinary incontinence
- Almost 50% of parous women suffer from some degree of pelvic organ prolapse

Conservative treatments such as pelvic floor muscle therapy (Kegel excercises) often fail because of patients' lack of compliance. Surgical options, although effective, suffer from a high rate of adverse effects and are typically a patient's last resort.

Revolutionary non-surgical laser treatment

Connective tissue in the vaginal walls is an important factor in pelvic organ support. Symptoms of pelvic floor dysfunction mainly arise from laxity in the vagina or its supporting ligaments, due to increasing age and vaginal childbirth.

Fotona's SMOOTH™ is a non-invasive non-ablative laser procedure for functional strengthening of connective tissue inside the vaginal wall, improving the pelvic floor support and diminishing symptoms of pelvic floor dysfunction.

What is SMOOTH™ mode?

SMOOTH™ gynecological procedures are based on the discovery that the delivery of an optimal sequence of heat pulses to the vaginal mucosa results in strengthening and rejuvenation of the vaginal wall. FotonaSmooth® is a gynecological laser system developed specifically for performing SMOOTH™ treatments.

FOTONA'S SMOOTH™ TREATMENTS INCLUDE:

IncontiLase® – Stress Urinary Incontinence

IntimaLase® – Vaginal Relaxation Syndrome

RenovaLase® – Vaginal Atrophy / Genitourinary Syndrome of Menopause

ProlapLase® – Pelvic Organ Prolapse

SIMPLE AND QUICK


- Ambulatory procedures
- 20-30 minutes/treatment

PATIENT FRIENDLY

- No anesthesia
- No downtime
- High patient satisfaction

Optimal Er:YAG laser wavelength

FotonaSmooth® operates at the optimal infrared Er:YAG laser wavelength (2.94 $\mu m)$ that coincides with the maximal absorption peak of the mucosal tissue. This ensures that the laser light is delivered to the tissue in a controlled superficial manner without the risk of affecting deeper-lying structures.




The Er:YAG wavelength sits at the highest peak of the absorbtion spectrum of water and is therefore completely absorbed within a few microns of mucosal tissue.

Unique SMOOTH™ mode


FotonaSmooth® delivers patented sequential Er:YAG SMOOTH™ mode laser pulses to the vaginal wall mucosa, generating controlled and optimal distribution of heat within the tissue, enabling collagen remodelling and neo-collagenesis.

The 2.94 µm wavelength, in conjunction with its patented SMOOTH™ mode delivery, allows for a highly controlled, safe procedure with no impact to any critical structures, including any penetration or disruption of the mucosal lining.

SMOOTH[™] mode pulse Optimal sequence of sub-ablative micro pulses



Unique sequential SMOOTH™ mode Er:YAG laser pulses generate an optimal structure of heat waves.



FotonaSmooth® – Non-surgical Er:YAG Laser for Vaginal and Pelvic Floor Health

Patient and physician-friendly

- Simple, fast in-office procedure
- Outstanding results, peer-reviewed
- No anesthesia required, non-ablative
- Minimal discomfort or downtime
- Excellent return on investment
- Multiple additional applications



Er:YAG treatments



Nd:YAG treatments

Easy to use

- Large touchscreen with intuitive graphical user interface
- Easily accessible presets for all applications
- Lightweight, ergonomic articulated-arm delivery system



Tools for gynecological treatments

G-SET™ - INTRAVAGINAL ACCESSORY SET



R11 full-beam titanium handpiece



PS03 patterned titanium handpiece



IncontiLase® and ProlapLase® treatment: 90° angular golden mirror titanium adapter



IntimaLase®, RenovaLase® and ProlapLase® treatment: 360° circular golden mirror titanium adapter



ocieai specululli



Wired laser speculum



FULL-BEAM AND PATTERNED HANDPIECES

- Zoom optics with spot sizes from 2-7 mm offers a wide range of treatments
- Collimated beam enables precise delivery of laser energy
- Titanium technology ensures robustness and durability
- Additional handpiece options for cervical treatments and soft tissue cutting





Wide Range of Treatments

FotonaSmooth® is more than a SMOOTH laser

Additional non-surgical Er:YAG gynecological procedures


The Er:YAG laser incorporated in FotonaSmooth® allows the user to perform other non-surgical treatments, such as:

- Vulvar depigmentation
- Labial and vulvar tightening

Surgical Er:YAG gynecological procedures

In additon, FotonaSmooth® can also be used for various surgical gynecological procedures, such as:

- Genital warts removal
- Cervical dysplasia vaporisation
- Labial trimming



Expand your practice

Additional Nd:YAG wavelength

FotonaSmooth® capabilities can be expanded with an additional Nd:YAG laser source. The Nd:YAG laser wavelength of 1064 µm exhibits the most homogeneous absorption in human tissue and thus represents a complementary wavelength to the Er:YAG wavelength with the shortest penetration depth in human tissue. Together, these golden-standard wavelengths are well suited for handling an exceptionally wide range of aesthetic and clinical procedures.

Key aesthetic treatments

Fr:YAG

Skin resurfacing, stretchmarks Scar revisions Lesions removal Acne scar treatments

Nd:YAG

Vascular lesion removal Removal of unwanted hair Wrinkle reduction Active acne treatments

Fotona Er:YAG and Nd:YAG laser system models

		Er:YAG GYNECOLOGY	Nd:YAG GYNECOLOGY	Er:YAG AESTHETICS	Nd:YAG AESTHETICS	
GYNECOLOGY MODELS						
FotonaSmooth® XS	20 W Er:YAG	YES	NO	Optional	NO	
FotonaSmooth® SP	20 W Er:YAG 30 W Nd:YAG	YES	YES	Optional	Optional	
FotonaSmooth® SP+	20 W Er:YAG 80 W Nd:YAG	YES	YES	Optional	Optional	
GENERAL USE MODELS						
XS Dynamis®	20 W Er:YAG	Optional	NO	YES	NO	
SP Spectro®	20 W Er:YAG 30 W Nd:YAG	Optional	Optional	YES	YES	
SP Dynamis®	20 W Er:YAG 80 W Nd:YAG	Optional	Optional	YES	YES	

FotonaSmooth® treatments



IncontiLase® A minimally invasive solution for stress urinary incontinence

How does IncontiLase® work?


- Fotona's 2.94 μm Er:YAG non-ablative laser with proprietary SMOOTH™ mode technology
- Improves urethral support by photothermal strengthening of the vaginal wall
- Works on connective tissue in the vaginal mucosa with emphasis on the anterior vaginal wall



Fotona SMOOTH™ mode treatment of the anterior vaginal wall

WHO IS ELIGIBLE FOR INCONTILASE®?


- The treatment works best in mild and moderate stress urinary incontinence patients, with very good results in severe stress urinary incontinence as well
- Mixed incontinence patients get relief in stress symptoms



Thermal camera images of the introitus show that SMOOTH™ pulses achieve the optimal peak temperature range (60-65°C) for collagen remodelling and initiation of neo-collagenesis.

Invasiveness Tissue	Non invasive methods	Invasive methods
Pelvic floor muscles	Kegel exercises Vaginal cones Electrostimulation	Surgical repair
Connective tissues	IncontiLase®	Anterior repair Sling procedures (TVT, TOT)

IncontiLase® treatment




One year after the treatment, significantly improved continence was reported in 77% of SUI patients.

Data presented in Lasers Surg Med. 47(9):689-97.

The effect of IncontiLase® therapy on the improvement of the grade of urinary incontinence (UI). Figure shows the distribution of patients (in %) with regard to the grade of incontinence (mild, moderate, severe, very severe) before treatment, at 2 months, 6 months and 1 year after the procedure.


Novel minimally invasive laser treatment of urinary incontinence in women. Ogrinc UB, Senčar S, Lenasi H. Lasers Surg Med., 2015: 47(9): 689-97.



Mild and moderate stress and mixed urinary incontinence



After IncontiLase® treatment



IncontiLase® treatment signifficantly reduced leakage quantity measured by 1h pad test.

Tien, Yi-Wen, et al. "Effects of laser procedure for female urodynamic stress incontinence on pad weight, urodynamics, and sexual function." International urogynecology journal 28.3 (2017): 469-476.



severe UI

One year following the IncontiLase® treatment, we found significant improvement in 77% of patients diagnosed with SUI.

Dr Sabina Sencar



Apart from non-invasiveness, the main advantage of IncontiLase® over surgery is that it can be applied as an ambulatory procedure, which means a lower economic burden.

Dr Urska Bizjak Ogrinc

Ogrinc UB, Senčar S, Lenasi H. Novel minimally invasive laser treatment of urinary incontinence in women. Lasers Surg Med. 47(9): 689-97.

FotonaSmooth® treatments



IntimaLase® A true incisionless laser treatment for vaginal relaxation syndrome


How does IntimaLase® work?

- Photothermally tightens the vaginal canal
- Mechanism of action is based on shrinking and thickening of the connective tissue in the vaginal wall

WHO IS ELIGIBLE FOR INTIMALASE®?

 Women with increased vaginal laxity due to childbirth and/or ageing

IntimaLase® treatment



Patients' assessment of sexual gratification improvement after IntimaLase® treatment.

Gaviria J, Lanz J. Laser Vaginal Tightening (LVT) – evaluation of a novel noninvasive laser treatment for vaginal relaxation syndrome. LAHA Journal of Laser and Health Academy, 2012(1); 46-58.

95% of patients reported an improvement of sexual gratification.

Patients' assessment of sexual gratification improvement after IntimaLase® treatment. Data presented in J. LA&HA, 2012(1); 46-58.



IntimaLase® treatment



95% of my patients assess vaginal tightness and sexual gratificatication as strongly or moderately improved after IntimaLase® treatment.

Dr Jorge Gaviria

Gaviria J, Lanz J. Laser Vaginal Tightening (LVT) – evaluation of a novel noninvasive laser treatment for vaginal relaxation syndrome. LAHA Journal of Laser and Health Academy, 2012(1); 46-58.




ProlapLase® Er:YAG SMOOTH™ mode laser treatment for pelvic organ prolapse

How does ProlapLase® work?

- Photothermal tightening of the tissue and contraction of the vaginal canal stimulating collagen remodeling and the synthesis of new collagen fibers
- A safe and non-invasive alternative to traditional methods
- Incisionless and virtually painless, with no cutting, bleeding or sutures
- Suitable also for higher grade prolapse

ProlapLase® treatment



The effect of ProlapLase® on cystocele grade distribution at baseline and follow-ups. The figure shows the effects of 2-7 laser treatments on 67 patients over a period of 3 years.

Erbium:YAG laser treatment of pelvic organ prolapses: 3 years follow-up. Sencar S., Bizjak-Ogrinc U., Vizintin

Significant reduction of average pelvic organ prolapse (POP) grade was already observed after the first treatment session. POP grade reduction continued to improve with additional treatments.

Data presented at IMS: 15th World Congress on Menopause.

WHO IS ELIGIBLE FOR PROLAPLASE®?

 Women suffering from pelvic organ prolapse



Pelvic organ prolapse



After ProlapLase® treatment



Dr Urska Bizjak Ogrinc



Dr Sabina Sencar

"In most of the patients the improvement occurred after the first treatment, with patients reporting better prolapse containment inside of the vaginal canal and less frequent occurrence of the prolapse falling out of vaginal canal."

Bizjak Obrinc U, Sencar S; Non-ablative laser treatment of cystocele It. J. Gynaecol. Obstet. 2017, 29: N. 120



RenovaLase® Gentle laser treatment of vaginal atrophy / genitourinary syndrome of menopause - GSM

How does RenovaLase® work?

- Non-ablative gentle photothermal treatment of the vaginal canal using very low energies that cause mild hyperthermia and induce microvascularisation and tissue regeneration without shrinking the collagen
- Restores normal vaginal mucosa structure and function
- Eliminates the need for long-term estrogen treatment

RenovaLase® treatment induces a significant improvement of genitourinary syndrome of menopause, including vaginal dryness and dyspareunia.
Additionaly, this treatment can be proposed in postmenopausal women who cannot be treated with hormones.

Dr Marco Gambacciani

Gambacciani M, Levancini M, Cervigni M. Vaginal erbium laser: the second-generation thermotherapy for the genitourinary syndrome of menopause. Climacteric. 2015, 18(5):757-763.




We observed significant improvement of symptoms and also regeneration of normal mucosal structure. Better and more lasting results compared to local estriol therapy.

Dr Adrian Gaspar


Gaspar A, Brandi H, Gomez V, Luque D. Efficacy of Erbium:YAG laser treatment compared to topical estriol treatment for symptoms of genitourinary syndrome of menopause. Lasers Surq Med. 2017, 49(2):160-168.

WHO IS ELIGIBLE FOR RENOVALASE®?

 Women suffering from vaginal atrophy symptoms such as vaginal dryness, irritation and dyspareunia




Atrophied vaginal mucosa




Vaginal mucosa after RenovaLase® treatment

Courtesy o

RenovaLase® treatment



RenovaLase® signifficantly decreased the severity of dyspareunia in breast cancer survivors suffering from GSM, the effect lasting more than 12 months after treatment.




Renovalase® had a significantly better effectiveness in relieving GSM symptoms.

RenovaLase® significantly improves GSM symptoms, including vaginal dryness and dyspareunia, in breast cancer survivors.

Gambacciani, Marco, and Marco Levancini. "Vaginal erbium laser as second-generation thermotherapy for the genitourinary syndrome of menopause: a pilot study in breast cancer survivors." Menopause 24.3 (2017): 316-319.

Comparison between RenovaLase® treatment and local estriol on severity of vaginal dryness up to 18 months after tretaments.


Gaspar, Adrian, et al. "Efficacy of Erbium: YAG laser treatment compared to topical estriol treatment for symptoms of genitourinary syndrome of menopause." Lasers in surgery and medicine 49.2 (2017): 160-168.



Non-ablative gentle photothermal treatment of the vaginal canal using mild hyperthermia via SMOOTH™ mode.



Enhanced synthesis of new proteoglycan and hialuronic acid molecules, which improves tissue water content and collagen structure.



The result is improved vascularisation and trophism of the vaginal wall and restored normal structure and function of the vaginal mucosa.

Scientifically proven results

Scientific research using FotonaSmooth®

STRESS URINARY INCONTINENCE

Effects of laser procedure for female urodynamic stress incontinence on pad weight, urodynamics, and sexual function. Tien YW, Hsiao SM, Lee CN, Lin HH. International Urogynecology Journal, 2017, 28(3):469-76.

Treatment of female stress urinary incontinence with Erbium-YAG laser in non-ablative mode. Pardo JI, Solà VR, Morales AA. Eur J Obstet Gynecol Reprod Biol. 2016 Sep;204:1-4.

Patient Satisfaction with Vaginal Erbium Laser Treatment of Stress Urinary Incontinence, Vaginal Relaxation Syndrome and Genitourinary Syndrome of Menopause. Barber MA, Eguiluz I. LAHA Journal of the Laser Health Academy, 2016(1):18-23.

Minimally invasive, non-ablative Er:YAG laser treatment of stress urinary incontinence in women – a pilot study. Fistonic N, Fistonic I, Findri Gustek S, Sorta Bilajac Turina I, Franic D, Vizintin Z, Kazic M, Hreljac I, Perhavec T, Lukac M. Lasers in Medical Science, 2016: vol. 31 (4) pp 635-43.

Novel minimally invasive laser treatment of urinary incontinence in women. Ogrinc UB, Senčar S, Lenasi H. Lasers Surg Med., 2015: 47(9): 689-97

Vaginal erbium laser: the second-generation thermotherapy for the genitourinary syndrome of menopause. Gambacciani M, Levancini M, Cerviqni M. Climacteric. 2015, 18(5):757-763.

Minimal Invasive Laser Treatment for Female Stress Urinary Incontinence. Khalafalla MM, Elbiaa AM, Abdelazim IA, Hussain M (2015). Obst&Gyn Int. Journal, 2015, 2(2) 00035.

First assessment of short-term efficacy of Er:YAG laser treatment on stress urinary incontinence in women: prospective cohort study. Fistonić N, Fistonić I, Lukanovič A, Guštek ŠF, Turina IS, Franić D. Climacteric. 2015 Oct;18 Suppl 1:37-42.

Minimally invasive laser procedure for early stages of stress urinary incontinence (SUI). Fistonić I, Findri-Guštek Š, Fistonić N. LAHA Journal of Laser and Health Academy, 2012(1); 67-74.

Novel Minimally Invasive VSP Er:YAG Laser Treatments in Gynecology. Vizintin Z, Rivera M, Fistonić I, Saracoğlu F, Guimares P, Gaviria J, Garcia V, Lukac M, Perhavec T, Marini L, LAHA Journal of Laser and Health Academy, 2012(1); 46-58.

VULVOVAGINAL ATROPHY / GENITOURINARY SYNDROME OF MENOPAUSE

Vaginal erbium laser as second-generation thermotherapy for the genitourinary syndrome of menopause: a pilot study in breast cancer survivors. Gambacciani M, Levancini M. Menopause, 2017, 24(3):316-319.

Efficacy of Erbium:YAG laser treatment compared to topical estriol treatment for symptoms of genitourinary syndrome of menopause. Gaspar A, Brandi H, Gomez V, Luque D. Lasers Surg Med. 2017, 49(2):160-168.

Treatment of Genitourinary Syndrome of Menopause with Erbium:YAG Laser: A Prospective Study of Efficacy and Safety of the Treatment for Women after Menopause of Natural Origin and Therapy-Induced Menopause in Breast Cancer Survivors. Bojanini JF. LAHA Journal of the Laser Health Academy, 2016(1):35-40.

Vaginal erbium laser: the second-generation thermotherapy for the genitourinary syndrome of menopause. Gambacciani M, Levancini M, Cervigni M. Climacteric. 2015, 18(5):757-763.

Vaginal Erbium Laser: the Second Generation Thermotherapy for the Genitourinary Syndrome of Menopause (GSM) in Breast Cancer Survivors. A preliminary report of a pilot study. Gambacciani M, Levancini M. It. J. Gynaecol. Obstet. 2015, 27: 9-11.

Short-term effect of vaginal erbium laser on the genitourinary syndrome of menopause. Gambacciani M, Levancini M. Minerva Ginecol. 2015 Apr;67(2):97-102.

Laser Treatment of Vaginal Atrophy in Post-menopause and Postgynecological Cancer Patients. Bojanini JF, Mejia AM. J. LA&HA, Vol. 2014, No.1; pp.65-71.

VAGINAL RELAXATION SYNDROME

Laser Vaginal Tightening with Non-ablative Er:YAG for Vaginal Relaxation Syndrome. Evaluation of Patient Satisfaction. Pardo JI, Solá Dalenz V. LAHA Journal of the Laser Health Academy, 2016(1):12-17.

Up to 3-Year Follow-up of Patients with Vaginal Relaxation Syndrome Participating in Laser Vaginal Tightening. Gaviria J, Korosec B, Fernandez J, Montero G. LAHA Journal of the Laser Health Academy, 2016(1):6-11.

Laser Vaginal Tightening (LVT) – evaluation of a novel noninvasive laser treatment for vaginal relaxation syndrome. Gaviria J, Lanz J. LAHA Journal of Laser and Health Academy, 2012(1); 46-58.

Novel Minimally Invasive VSP Er:YAG Laser Treatments in Gynecology. Vizintin Z, Rivera M, Fistonić I, Saracoğlu F, Guimares P, Gaviria J, Garcia V, Lukac M, Perhavec T, Marini L, LAHA Journal of Laser and Health Academy, 2012(1); 46-58.

PELVIC ORGAN PROLAPSE

Non-ablative vaginal erbium YAG laser for the treatment of cystocele. Bizjak Ogrinc U, Sencar S. It. J. Gynaecol. Obstet. 2017, 29:19-25.

Laser Treatment of Higher Grade Cystocele. Bizjak Ogrinc U, Sencar S. LAHA Journal of the Laser and Health Academy, 2013 (1); S22.

Erbium Laser in Gynecology. Vizintin Z, Lukac M, Kazic M, Tettamanti M. Climacteric. 2015, 18(1):4-8.

Training & Education

Regular clinical workshops – make a SMOOTH start!

Comprehensive workshops

Training is provided in cooperation with the Laser and Health Academy under the guidance of experts in medical laser technology.

The extensive workshops, where participants engage in live demonstrations and gain in-depth uderstanding of laser physics and laser tissue interaction, provide the needed insight into the fundamentals of non-invasive gynecological treatments and other procedures that can be perfomed with the FotonaSmooth® system.



Fotona and the Laser & Health Academy

Fotona has partnered with the Laser & Health Academy (LA&HA) to help support the professional growth of medical practitioners. To get the most out of your Fotona laser system, our practitioner workshops, co-organized with LA&HA (www.laserandhealth.com), provide hands-on demonstrations of our lasers by international clinical experts.



ANNUAL INTERNATIONAL LASER & HEALTH ACADEMY SYMPOSIUM

- Attended by several hundred physicians every year
- Newest research and treatments in different fields of laser medicine, including gynecology
- For more information, contact info@laserandhealth.com

THE WORKSHOPS COVER:

- Laser safety and physics
- Laser-tissue interaction
- Extensive theoretical and hands-on application training
- A visit to a clinic for live patient demonstrations





Committed to Engineering

The Highest Performance, Best Made Laser Systems in the World

since 1964

Founded in 1964, only four years after the invention of the very first laser, Fotona is one of the most experienced developers of high-technology laser systems. Fotona today is a world-leading medical laser company recognized for its innovative, award-winning laser systems for applications in gynecology, dentistry, surgery and aesthetics & dermatology. Based in the US and EU, with corporate headquarters in Dallas, Texas, and Ljubljana, Slovenia, Fotona's business philosophy is to continuously choose perfection to ensure the maximum performance and efficacy of its medical devices.



Fotona, d. o. o. Stegne 7 1000 Ljubljana Slovenia EU Fotona, LLC 2307 Springlake Road #518 Dallas, TX 75234 USA Fotona GmbH Hohlbachweg 2 73344 Gruibingen, Germany EU Fotona Beauty Light, (Suzhou) Medical Devices Co, Ltd. No 2, Zengfu Road, Guli Town Changshu City, Jiangsu Province CHINA, 215515

info@fotona.com www.fotona.com

All Fotona medical lasers are CE marked and approved to be sold in the EU. For countries where specific national approvals or clearances are required, some of the products and/or applications may not yet have been approved. Please check with Fotona, your local Fotona distributor or your national regulatory body about whether a specific product or application has been approved to be marketed and sold in your country.



